

EXCELLENT PERFORMANCE AND SAFETY!

Pre & post vacuum tabletop “Class B” autoclave, available in two sizes; 18 litres and 23 litres.

The autoclaves require minimal bench space making it the ideal sterilisation solution for veterinary clinics.



SAFETY & MONITORING

Double safety protection system prevents the cycle from starting if the door is not correctly locked. This system also prevents the door from opening if the pressure inside the chamber is not equal to atmospheric pressure outside the chamber.

Pressure safety valve prevents over pressurisation in the chamber and steam generator.

Automatic power cut off if a short circuit or electrical fault occurs.

Able to detect and identify the exact cause of any problem, and to issue a specific error code to the operator.



EXCELLENT PERFORMANCE

- Full automatic table top autoclave with pre-set programs.
- Complies with European standard EN13060.
- Work very quietly.
- In built independent rapid steam generator ensures quick sterilisation cycles.
- Rapid post sterilisation vacuum drying.
- With one clean water tank and one used water tank.
- LCD screen displays time, temperature, pressure, process alerts and instrument conditions.
- User friendly interface.
- Water level switch in the main tank controls the min and max water levels.
- Automatic warning for maintenance.



	ICANCLAVE D218	ICANCLAVE D223
CenQuip Code	STE-18-D	STE-23-D
Chamber Dimension	Ø247x350 mm	Ø247x450 mm
Chamber Volume	18 Litres	23 Litres
Number of Trays	3 Standard. Max. 4	3 Standard. Max. 4
Material of Trays	Aluminium	Aluminium
Size of Trays.	192x285mm	192x380mm
Voltage Freq.	220/110V - 50/60Hz	220/110V - 50/60Hz
POWER.r	2000W	2200W
Overall Dimension (WxHxD)	480 x 450 x 600 mm	480 x 455 x 700 mm
Weight	47 kgs	51 kgs