



**CenQuip**<sup>®</sup>



# EXAPad

EXCELLENCE AT YOUR FINGERTIPS



**ECHOTouch**



**ECHOVoice**



**ECHOPad**

For more information  
Phone: 02 9679 5733  
Email: [info@cenquip.com.au](mailto:info@cenquip.com.au)

 **ECM**  
Echo Control Medical

## A new way of approaching ultrasound scanning.

The streamlined interface provides every user with a simple and completely intuitive use of the ultrasound scanner. The new applications, EchoTouch and EchoVoice optimise the efficiency making it a handy ultrasound scanner.

A high resolution 15" touchscreen, with unequalled robustness for quality imaging, battery life of more than two hours, exceptional lightweight and handling make the EXAPAD a high end portable scanner capable of answering the demands veterinarians.

Vet specific algorithms make diagnostics easier. Patient files, reports, annotations, specific measurements, are all adapted to veterinary practice.



**ECHOTouch:** Thanks to the new EchoTouch interface system, only the necessary functions are present on the screen, making use faster and more efficient. The advanced settings remain for simple access as well.



**ECHOVoice:** The embedded software, EchoVoice, gives you a voice control over the main functions of the scanner, allowing you to work hands free, even when you are further away from the device.



**ECHOPad:** Remote access with EchoPad installed on your Ipad, you have access to a second command board for your Exapad, allowing you to pilot your ultrasound scanner remotely. You can save your own presets in your iPad that can be used in the Exapad.

### Technical Specifications

- ✓ **Screen:** 15" Touchscreen
- ✓ **Weight:** 6.5kg
- ✓ **Autonomy:** > 2 hours
- ✓ Shockproof transport case
- ✓ Color Doppler CFM
- ✓ **Imaging modes:**
  - Mode B, M, B/B, B/M
- ✓ Pulse wave doppler PW
- ✓ Wi-Fi, Bluetooth
- ✓ Kensington security



• Sales • Service • Installation & Training  
• Technical Support

For more information

Phone: 02 9679 5733 or Email: [info@cenquip.com.au](mailto:info@cenquip.com.au)